

---

## Product Datasheet

### TLR8 Antibody (orb1238754)

**Description**

TLR8 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested**

ELISA, ICC, IF, WB

**Applications**
**Immunogen**

TLR8 antibody was raised against a peptide corresponding to 16 amino acids near the middle of human TLR8. The immunogen is located within amino acids 790 - 840 of TLR8.

**Target**

TLR8

**Preservatives**

TLR8 Antibody is supplied in PBS containing 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

TLR8 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

TLR8 antibody can be used for detection of TLR8 by Western blot at 0.5 to 2 µg/mL. Antibody can also be used for immunocytochemistry starting at 2 µg/mL. For immunofluorescence start at 10 µg/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

**Isotype**

IgG

**Clonality**


Polyclonal

**MW**

Predicted: 112 kDa Observed: 107 kDa

**Uniprot ID**
[Q9NR97](#)
**NCBI**
[NP\\_619542](#)
**Dilution Range**


TLR8 antibody can be used for detection of TLR8 by Western blot at 0.5 to 2 µg/mL. Antibody can also be used for immunocytochemistry starting at 2 µg/mL. For immunofluorescence start at 10 µg/mL. Antibody validated: Western Blot in human samples;



Western blot analysis of TLR8 in Daudi c...



Immunocytochemistry of TLR8 in Daudi cel...



Immunofluorescence of TLR8 in Daudi cell...