
Product Datasheet

Tlr11 Antibody (orb1238766)

Description

Tlr11 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, WB

Applications
Immunogen

TLR11 antibody was raised against a peptide corresponding to 16 amino acids between the middle and carboxy terminus of mouse TLR11. The immunogen is located within the last 50 amino acids of TLR11.

Target

Tlr11

Preservatives

TLR11 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TLR11 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

TLR11 antibody polyclonal antibody can be used for detection of TLR11 by Western blot at 0.5 to 2 µg/mL. Antibody can also be used for immunocytochemistry starting at 10 µg/mL. For immunofluorescence start at 10 µg/mL. Antibody validated: Western Blot in mouse samples; Immunocytochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.

Isotype


IgG

Clonality


Polyclonal

Uniprot ID
[Q6R5P0](#)
NCBI
[NP_991388](#)
Dilution Range


TLR11 antibody polyclonal antibody can be used for detection of TLR11 by Western blot at 0.5 to 2 µg/mL. Antibody can also be used for immunocytochemistry starting at 10 µg/mL. For immunofluorescence start at 10 µg/mL. Antibody validated: Western Blot in mouse



Western blot analysis of TLR11 in RAW264...



Immunocytochemistry of TLR11 in RAW264.7...



Immunofluorescence of TLR11 in Raw264 ce...