

Product Datasheet

TIRAP Antibody (orb1238771)

Description

TIRAP Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

TIRAP antibody was raised against a peptide corresponding to 15 amino acids near the C-terminus of mouse TIRAP. The immunogen is located within the last 50 amino acids of TIRAP.

Target

TIRAP

Preservatives

TIRAP Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TIRAP antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

TIRAP antibody can be used for detection of TIRAP by Western blot at 1 to 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 2 µg/mL. For immunofluorescence start at 10 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 24, 26 kDa Observed: 24 kDa

Uniprot ID


[P58753](#)

NCBI


[NP_473437](#)

Dilution Range


TIRAP antibody can be used for detection of TIRAP by Western blot at 1 to 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 2 µg/mL. For immunofluorescence start at 10 µg/mL. Antibody validated: Western Blot in human samples;



Western blot analysis of TIRAP in human ...



Immunohistochemical staining of human he...



Immunofluorescence of TIRAP in Human Hea...