

Product Datasheet

CREB1 antibody (orb12388)

Description Rabbit polyclonal antibody to CREB1

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Recombinant fusion protein of human CREB1(NP_604391.1).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 43kDa


Uniprot ID **P16220**

Dilution Range WB: 1:500-1:2000, IF/ICC: 1:50-1:100


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of mouse kidney



Immunohistochemical staining of mouse brain



Immunohistochemical staining of human lung