



Product Datasheet

StreptII-tag Antibody (orb1238868)

Description	StreptII-tag Antibody
Species/Host	Rabbit
Reactivity	Tag
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Anti-StreptII-tag antibody (orb1238868) StreptII-tag antibody was raised against the amino acid sequence "asWSHPQFEKga".
Preservatives	StreptII-tag Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	StreptII-tag antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB: 0.25-1 µg/mL Antibody validated: StreptII-tag antibody can be used for the detection of StreptII-tagged proteins
Expiration Date	12 months from date of receipt.



Western Blot Validation with StreptII-tag...



Western Blot Validation with StreptII-tag...