
Product Datasheet

SNRPN Antibody (orb1238918)

Description

SNRPN Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

SNRPN antibody was raised against a 20 amino acid peptide near the amino terminus of human SNRPN. The immunogen is located within amino acids 30 - 80 of SNRPN.

Target

SNRPN

Preservatives

SNRPN antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SNRPN antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

SNRPN antibody can be used for detection of SNRPN by Western blot at 0.5 - 1 µg/mL. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality


Polyclonal

MW

Predicted: 26 kDa Observed: 25 kDa

Uniprot ID
[P63164](#)
NCBI
[NP_003088](#)
Dilution Range

SNRPN antibody can be used for detection of SNRPN by Western blot at 0.5 - 1 µg/mL. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.




Western blot analysis of SNRPN in THP-1 ...



Immunohistochemistry of SNRPN in human

S...



Immunofluorescence of SNRPN in human spl...