

Product Datasheet

EBI3 Antibody (orb1238953)

Description

EBI3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

Anti-EBI-3 antibody (orb1238953) was raised against a peptide corresponding to 16 amino acids near the center of mature human EBI-3. The immunogen is located within amino acids 150 - 200 of EBI3.

Target

EBI3

Preservatives

EBI3 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

EBI3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 25kD
Observed: 30 kD (25kD + 2 N-Glycosylations)

Uniprot ID

[Q14213](#)

NCBI


[NP_005746](#)

Dilution Range


WB: 0.5 - 2 µg/mL; IF: 20 µg/mL ; IHC: 2.5 µg/mL;. Antibody validated: Western Blot in human and rat samples; Immunofluorescence and Immunohistochemistry in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot
Validation in Rat
Spleen Ti...



Recombinant Protein
Validation. Loading:...



Immunohistochemistry
Validation of EBI-3...