
Product Datasheet

Non-structural protein 9 Antibody (orb1239069)

Description

Non-structural protein 9 Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Immunogen

Anti-SARS-CoV-2 (COVID-19) NSP9 antibody (orb1239069) was raised against a peptide corresponding to 13 amino acids near the center of SARS-CoV-2 (COVID-19) NSP9 protein. The immunogen is located within 30-80 amino acids of the SARS-CoV-2 (COVID-19) NSP9 protein.

Target

Non-structural protein 9

Preservatives

SARS-CoV-2 (COVID-19) NSP9 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 (COVID-19) NSP9 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality


Polyclonal

Uniprot ID
[P0DTC1](#)
NCBI
[YP_009742616.1](#)
Dilution Range


IHC: 0.5 µg/ml. Antibody validated: Immunohistochemistry in human samples. SARS-CoV-2 (COVID-19) NSP9 antibody can detect 2 ng of free peptide at 1 µg/mL in ELISA. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Immunohistochemistry
Validation of SARS-...



ELISA Validation.
Antibodies: SARS-CoV-2...