


Product Datasheet

HEL Antibody (orb1239078)


Description	HEL Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Anti-SARS-CoV-2 (COVID-19) NSP13 (Helicase)antibody (orb1239078) was raised against a peptide corresponding to 14 amino acids near the center of the SARS-CoV-2 (COVID-19) NSP13 (Helicase) protein. The immunogen is located between 440-490 amino acids of SARS-CoV-2 (COVID-19) NSP13 (Helicase).
Target	HEL
Preservatives	SARS-CoV-2 (COVID-19) NSP13 (Helicase) Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	SARS-CoV-2 COVID-19) NSP13 (Helicase) antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	PODTC1
NCBI	YP_009725308.1
Dilution Range	IHC: 0.5 µg/mL. Antibody validated: Immunohistochemistry in human samples. SARS-CoV-2 (COVID-19) NSP13 antibody can detect 2 ng of free peptide at 1 µg/mL in ELISA. It can detect SARS-CoV-2 NSP13 recombinant protein by ELISA. All other applications and species not yet tested
Expiration Date	12 months from date of receipt.



Immunohistochemistry Validation of SARS-...



ELISA Validation with SARS-CoV-2 (COVID-...



ELISA Validation. Antibodies: SARS-CoV-2...