

---

## Product Datasheet

### PIAS4 Antibody (orb1239234)

**Description**

PIAS4 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, IHC-P, WB

**Immunogen**

PIAS4 antibody was raised against a 19 amino acid synthetic peptide near the amino terminus of human PIAS4. The immunogen is located within amino acids 50 - 100 of PIAS4.

**Target**

PIAS4

**Preservatives**

PIAS4 Antibody is supplied in PBS containing 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

PIAS4 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

PIAS4 antibody can be used for detection of PIAS4 by Western blot at 1 µg/mL. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in rat samples. All other applications and species not yet tested.

**Isotype**

IgG

**Clonality**


Polyclonal

**MW**


Predicted: 56 kDa Observed: 56 kDa

**Uniprot ID**
[Q8N2W9](#)
**NCBI**
[NP\\_056981](#)
**Dilution Range**

PIAS4 antibody can be used for detection of PIAS4 by Western blot at 1 µg/mL. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in rat samples. All other



Western blot analysis of PIAS4 in HL60 c...



Immunohistochemistry of PIAS4 in rat bra...