

Product Datasheet

ANP32E Antibody (orb1239242)

Description

ANP32E Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

PHAP III antibody was raised with a synthetic peptide corresponding to amino acids close to carboxy terminus of human PHAP III. The immunogen is located within amino acids 200 - 250 of PHAP III.

Target

ANP32E

Preservatives

PHAP III Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PHAP III antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures..

Note

For research use only

Application notes

PHAP III antibody can be used for detection of PHAP III by Western blot at 1 µg/mL. A band at approximately 35 kDa can be detected. Antibody can also be used for immunohistochemistry starting at 2 µg/mL. For immunofluorescence start at 10 µg/mL. Antibody validated: Western Blot in human and rat samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

35 kDa

Uniprot ID


[Q9BTT0](#)

NCBI


[NP_112182](#)

Dilution Range


PHAP III antibody can be used for detection of PHAP III by Western blot at 1 µg/mL. A band at approximately 35 kDa can be detected. Antibody can also be used for immunohistochemistry starting at 2 µg/mL. For immunofluorescence start at 10 µg/mL.



Western blot analysis of PHAP III expres...



Immunohistochemistry of PHAP III in huma...



Immunofluorescence of PHAP III in Human

...