

---

## Product Datasheet

### NIPSNAP3B Antibody (orb1239365)

**Description**

NIPSNAP3B Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, IF, IHC-P, WB

**Immunogen**

NIPSNAP3B antibody was raised against a 16 amino acid synthetic peptide near the center of human NIPSNAP3B. The immunogen is located within amino acids 40 - 90 of NIPSNAP3B.

**Target**

NIPSNAP3B

**Preservatives**

NIPSNAP3B Antibody is supplied in PBS containing 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

NIPSNAP3B antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

NIPSNAP3B antibody can be used for detection of NIPSNAP3B by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 2.5 µg/mL. For immunofluorescence start at 5 µg/mL. Antibody validated: Western Blot in mouse samples; Immunohistochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.

**Isotype**


IgG

**Clonality**


Polyclonal

**Uniprot ID**
[Q9BS92](#)
**NCBI**
[Q9BS92](#)
**Dilution Range**


NIPSNAP3B antibody can be used for detection of NIPSNAP3B by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 2.5 µg/mL. For immunofluorescence start at 5 µg/mL. Antibody



Western blot analysis of NIPSNAP3B in mo...



Immunohistochemistry of NIPSNAP3B in mou...



Immunofluorescence of NIPSNAP3B in mouse...