
Product Datasheet

MYH3 Antibody (orb1239414)

Description

MYH3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IHC-P, WB

Immunogen

MYH3 antibody was raised against a 19 amino acid peptide near the amino terminus of human MYH3. The immunogen is located within amino acids 30 - 80 of MYH3.

Target

MYH3

Preservatives

MYH3 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

MYH3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

MYH3 antibody can be used for detection of MYH3 by Western blot at 1 µg/ml. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 240 kDa Observed: 245 kDa

Uniprot ID
[P11055](#)
NCBI
[NP_002461](#)
Dilution Range


MYH3 antibody can be used for detection of MYH3 by Western blot at 1 µg/ml. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in mouse samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of MYH3 in fetal h...



Immunohistochemistry of MYH3 in mouse sk...