


Product Datasheet

IL17A Antibody (orb1239477)


Description	IL17A Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	Anti-IL-17 antibody (orb1239477) was raised against a peptide corresponding to 16 amino acids near the center of mature human IL-17. The immunogen is located within amino acids 50-100 of IL-17.
Target	IL17A
Preservatives	IL-17 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	IL-17 antibody can be stored at 4°C for up to one year. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 17.5kDa Observed: 18 kDa
Uniprot ID	Q16552
NCBI	NP_002181
Dilution Range	WB: 0.125-2 µg/mL; IF: 5-20 µg/mL; IHC: 1-5 µg/mL .Antibody validated: Western Blot in human samples; Immunofluorescence in mouse samples; Immunohistochemistry in human and mouse samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western Blot Validation with Recombinant...



Western Blot Validation with Human Splee...



Immunohistochemistry Validation of IL-17...