

Product Datasheet

KDM6B Antibody (orb1239526)

Description

KDM6B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, ICC, IF, WB

Immunogen

Anti-JMJD3 antibody (orb1239526) was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of human JMJD3. The immunogen is located within amino acids 1370-1420 of JMJD3.

Target

KDM6B

Preservatives

JMJD3 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

JMJD3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 180Observed: 180 kD

Uniprot ID

[O15054](#)

NCBI


[EAW90126](#)

Dilution Range


WB: 0.25-4 µg/mL; ICC: 2.5 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunocytochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot Validation in Human, Mouse ...



Independent Antibody Validation (IAV) vi...



Western Blot Validation in HeLa Cell Lys...