


Product Datasheet

IL32 Antibody (orb1239659)

Description	IL32 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Anti-IL-32 antibody (orb1239659) was raised against a peptide corresponding to 16 amino acids near the center of human IL-32. The immunogen is located within amino acids 50 - 100 of IL-32.
Target	IL32
Preservatives	IL-32 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	IL-32 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 26kD Observed: 26 kD
Uniprot ID	P24001
NCBI	AAH09401
Dilution Range	WB: 0.5-5 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western Blot Validation with Recombinant...



Western Blot Validation in Jurkat Cell L...