

---

## Product Datasheet

### IFIT3 Antibody (orb1239695)

**Description**

IFIT3 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, IHC-P, WB

**Immunogen**

IFIT3 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human IFIT3. The immunogen is located within amino acids 380 - 430 of IFIT3 .

**Target**

IFIT3

**Preservatives**

IFIT3 Antibody is supplied in PBS containing 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

IFIT3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

**Note**

For research use only

**Application notes**

IFIT3 antibody can be used for detection of IFIT3 by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry at 5 µg/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in human samples. All other applications and species not yet tested.

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**


Predicted: 54 kDa Observed: 50 kDa

**Uniprot ID**
[O14879](#)
**NCBI**
[NP\\_001026853](#)
**Dilution Range**


IFIT3 antibody can be used for detection of IFIT3 by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry at 5 µg/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in human samples. All other applications and species not yet tested.

**Expiration Date**

12 months from date of receipt.



Western blot analysis of IFIT3 in HepG2 ...



Immunohistochemistry of IFIT3 in human l...