


---


## Product Datasheet

### HOOK3 Antibody (orb1239706)


<b>Description</b>	HOOK3 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC-P, WB
<b>Immunogen</b>	HOOK3 antibody was raised against a 19 amino acid peptide from near the amino terminus of human HOOK3. The immunogen is located within amino acids 210 - 260 of HOOK3.
<b>Target</b>	HOOK3
<b>Preservatives</b>	HOOK3 Antibody is supplied in PBS containing 0.02% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	HOOK3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
<b>Note</b>	For research use only
<b>Application notes</b>	HOOK3 antibody can be used for the detection of HOOK3 by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry at 10 µg/ml. Antibody validated: Western Blot in mouse samples; Immunohistochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	Predicted: 79 kDa Observed: 77 kDa
<b>Uniprot ID</b>	<a href="#">Q86VS8</a>
<b>NCBI</b>	<a href="#">NP_115786</a>
<b>Dilution Range</b>	HOOK3 antibody can be used for the detection of HOOK3 by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry at 10 µg/ml. Antibody validated: Western Blot in mouse samples; Immunohistochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HOOK3 in mouse ...



Immunohistochemistry of HOOK3 in mouse b...



Immunofluorescence of HOOK3 in mouse bra...