



Product Datasheet

GAG Antibody (orb1239710)


Description	GAG Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Mouse monoclonal HIV-1 p24 antibody was raised against a recombinant full-length HIV-1 p24 protein.
Target	GAG
Preservatives	HIV-1 p24 Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	HIV-1 p24 Monoclonal Antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	orb1239710 can be used for detection of p24 by Western blot or ELISA at 0.2 - 0.5µg/mL.
Isotype	IgG1
Clonality	Monoclonal
MW	Observed: 24 kDa, 41 kDa, and 55 kDa
Uniprot ID	P04591
NCBI	AAB50258
Dilution Range	orb1239710 can be used for detection of p24 by Western blot or ELISA at 0.2 - 0.5µg/mL.
Expiration Date	12 months from date of receipt.



Western blot analysis of 20 ng of (A) vi...



Titration ELISA of anti-HIV-1 p24 antibo...



Sandwich ELISA analysis of recombinant p...