




Product Datasheet

FAF1 Antibody (orb1239816)

Description	FAF1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IHC-P, WB	
Immunogen	FAF1 antibody was raised against a 19 amino acid synthetic peptide from near the carboxy terminus of human FAF1. The immunogen is located within the last 50 amino acids of FAF1.	 Western blot analysis of FAF1 in THP-1 c...
Target	FAF1	 Immunohistochemistry of FAF1 in rat sple...
Preservatives	FAF1 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	FAF1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	 Immunofluorescence of FAF1 in rat spleen...
Note	For research use only	
Application notes	FAF1 antibody can be used for detection of FAF1 by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 2.5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 74 kDa Observed: 80 kDa	
Uniprot ID	Q9UNN5	
NCBI	NP_008982	
Dilution Range	FAF1 antibody can be used for detection of FAF1 by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 2.5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human	