
Product Datasheet

STK17A Antibody (orb1239863)

Description

STK17A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, WB

Applications
Immunogen

DRAK1 antibody was raised against a peptide corresponding to amino acids near the amino terminus of human DRAK1. The immunogen is located within the first 50 amino acids of DRAK1.

Target

STK17A

Preservatives

DRAK1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

DRAK1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

DRAK1 antibody can be used for detection of DRAK1 by Western blot at 1 µg/mL. An approximately 50 kDa band can be detected. Antibody can also be used for immunocytochemistry starting at 2 µg/mL. For immunofluorescence start at 10 µg/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality


Polyclonal

MW


50 kDa

Uniprot ID
[Q9UEE5](#)
NCBI
[Q9UEE5](#)
Dilution Range


DRAK1 antibody can be used for detection of DRAK1 by Western blot at 1 µg/mL. An approximately 50 kDa band can be detected. Antibody can also be used for immunocytochemistry starting at 2 µg/mL.



Western blot analysis of DRAK1 in (A) MO...



Immunocytochemistry staining of MOLT4 ce...



Immunofluorescence of DRAK1 in Molt cell...