

Product Datasheet

DDX41 Antibody (orb1239893)

Description

DDX41 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

DDX41 antibody was raised against a 17 amino acid peptide near the amino terminus of human DDX41. The immunogen is located within amino acids 230 - 280 of DDX41.

Target

DDX41

Preservatives

DDX41 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

DDX41 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

DDX41 antibody can be used for detection of DDX41 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 68 kDa Observed: 70 kDa

Uniprot ID

[Q9UJV9](#)

NCBI


[NP_057306](#)

Dilution Range


DDX41 antibody can be used for detection of DDX41 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western blot analysis of DDX41 in rat br...



Immunohistochemistry of DDX41 in rat bra...



Immunofluorescence of DDX41 in rat brain...