

Product Datasheet

TET1 Antibody (orb1240105)

Description

TET1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

Anti-TET1 antibody (orb1240105) was raised against an 18 amino acid peptide near the carboxy terminus of human TET1. The immunogen is located within amino acids 2030 - 2080 of TET1.

Target

TET1

Preservatives

TET1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TET1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 235 kDa Observed: 235 kDa

Uniprot ID

[Q8NFU7](#)

NCBI


[NP_085128](#)

Dilution Range


WB: 1 µg/mL; IHC: 1 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot and Immunohistochemistry in human, mouse and rat samples; Immunocytochemistry in human samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



WB Validation of TET1 in Human Cell Line...



WB Validation of TET1 in Mouse Tissues. ...



WB Validation of TET1 in Rat Tissues. Lo...