
Product Datasheet

ARHGAP18 Antibody (orb1240207)

Description

ARHGAP18 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, WB

Applications
Immunogen

ARHGAP18 antibody was raised against a 15 amino acid synthetic peptide near the carboxy terminus of human ARHGAP18. The immunogen is located within the last 50 amino acids of ARHGAP18.

Target

ARHGAP18

Preservatives

ARHGAP18 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

ARHGAP18 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

ARHGAP18 antibody can be used for detection of ARHGAP18 by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunocytochemistry starting at 2.5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in mouse samples; Immunocytochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.

Isotype

IgG


Clonality

Polyclonal


Uniprot ID
[Q8N392](#)
NCBI
[NP_277050](#)
Dilution Range

ARHGAP18 antibody can be used for detection of ARHGAP18 by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunocytochemistry starting at 2.5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in mouse samples; Immunocytochemistry in mouse samples

A B




Western blot analysis of ARHGAP18 in 3T3...



Immunocytochemistry of ARHGAP18 in 3T3

C...



Immunofluorescence of ARHGAP18 in 3T3

ce...