


Product Datasheet

TNFSF12 Antibody (orb1240271)


Description	TNFSF12 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	Anti-TWEAK antibody (orb1240271) was raised against recombinant human TWEAK protein.
Target	TNFSF12
Preservatives	TWEAK Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	TWEAK antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 27 and 36 kDa Observed: 27 and 36 kD
Uniprot ID	O43508
NCBI	NP_003800
Dilution Range	WB: 0.5-4 µg/mL; IHC: 5-10 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunohistochemistry in human and mouse samples; Immunofluorescence in mouse samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western Blot Validation in HeLa Cell Lys...



Western Blot Validation with Recombinant...



Western Blot Validation in Mouse Pancrea...