


Product Datasheet

VCL Antibody (orb1240274)


Description	VCL Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Anti-Vinculin antibody (orb1240274) was raised against a peptide corresponding to 23 amino acids near the center of human Vinculin.
Target	VCL
Preservatives	Vinculin antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Vinculin antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 124kD Observed: 124kD
Uniprot ID	P18206
NCBI	NP_003364
Dilution Range	WB: 0.5-2 µg/mL; ICC: 5 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunocytochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Independent Antibody Validation (IAV) vi...



Western Blot Validation in Human, Mouse ...



Western Blot Validation in PC3 Cell Lysa...