


## Product Datasheet

# Bebtelovimab (LY-CoV1404/LY3853113) Antibody, Monoclonal (orb1240314)

<b>Description</b>	Bebtelovimab (LY-CoV1404/LY3853113) Antibody, Monoclonal
<b>Species/Host</b>	Human
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1-kappa
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.



SDS-PAGE image