
Product Datasheet

H1FOO Antibody (orb1240387)

Description	H1FOO Antibody	
Species/Host	Rabbit	
Reactivity	Human	40kDa
Conjugation	Unconjugated	30kDa
Tested Applications	ELISA, WB	20kDa
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human H1FOO.	Antibody used in WB on Human HeLa at 0.2...
Target	H1FOO	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store H1FOO antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	H1FOO antibody can be used for detection of H1FOO by ELISA at 1:312500. H1FOO antibody can be used for detection of H1FOO by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	36 kDa	
Uniprot ID	Q81ZA3	
NCBI	NP_722575	
Dilution Range	H1FOO antibody can be used for detection of H1FOO by ELISA at 1:312500. H1FOO antibody can be used for detection of H1FOO by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	