
Product Datasheet

TCF12 Antibody (orb1240511)

Description	TCF12 Antibody	910 _u
Species/Host	Rabbit	910 _u
Reactivity	Human, Mouse, Rat	910 _u
Conjugation	Unconjugated	910 _u
Tested Applications	ELISA, WB	910 _u
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TCF12.	Antibody used in WB on Human Lung at 0.0...
Target	TCF12	910 _u
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	910 _u
Form/Appearance	Liquid	910 _u
Concentration	batch dependent	910 _u
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TCF12 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Hum. Fetal Brain ...
Note	For research use only	910 _u
Application notes	TCF12 antibody can be used for detection of TCF12 by ELISA at 1:62500. TCF12 antibody can be used for detection of TCF12 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	910 _u
Clonality	Polyclonal	910 _u
MW	73 kDa	910 _u
Uniprot ID	Q99081	910 _u
NCBI	NP_003196	910 _u
Dilution Range	TCF12 antibody can be used for detection of TCF12 by ELISA at 1:62500. TCF12 antibody can be used for detection of TCF12 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	910 _u
Expiration Date	12 months from date of receipt.	910 _u

910_u910_u910_u910_u

Antibody used in WB on Human Lung at 0.0...

TCF12

90

65

40

31

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22

22