


Product Datasheet

SREBF2 Antibody (orb1240596)

Description	SREBF2 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on mouse glomerular ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SREBF2.	
Target	SREBF2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SREBF2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Human Raji cells ...
Note	For research use only	
Application notes	SREBF2 antibody can be used for detection of SREBF2 by ELISA at 1:62500. SREBF2 antibody can be used for detection of SREBF2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	124 kDa	
Uniprot ID	Q12772	
NCBI	NP_004590	
Dilution Range	SREBF2 antibody can be used for detection of SREBF2 by ELISA at 1:62500. SREBF2 antibody can be used for detection of SREBF2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	