

Product Datasheet

HIST3H3 Antibody (orb1240657)

Description	HIST3H3 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide directed towards the N terminal of human HIST3H3	
Target	HIST3H3	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store H3F3A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Human Brain at 1 ...
Note	For research use only	
Application notes	H3F3A antibody can be used for detection of H3F3A by ELISA at > 1:3000. H3F3A antibody can be used for detection of H3F3A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	15 kDa, 15 kDa, 16 kDa, 15 kDa, 23 kDa	
Uniprot ID	P06351	
NCBI	NP_002098	
Dilution Range	H3F3A antibody can be used for detection of H3F3A by ELISA at > 1:3000. H3F3A antibody can be used for detection of H3F3A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	