

Product Datasheet

HUS1B Antibody (orb1240766)

Description	HUS1B Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Liver at 0...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HUS1B.	
Target	HUS1B	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store HUS1B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	HUS1B antibody can be used for detection of HUS1B by ELISA at 1:312500. HUS1B antibody can be used for detection of HUS1B by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	31 kDa	
Uniprot ID	Q8NHY5	
NCBI	NP_683762	
Dilution Range	HUS1B antibody can be used for detection of HUS1B by ELISA at 1:312500. HUS1B antibody can be used for detection of HUS1B by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	