
Product Datasheet

BMPER Antibody (orb1240767)

Description

BMPER Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human BMPER.

Antibody used in WB on Human Fetal heart...

Target

BMPER

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store BMPER antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in IHC on Human Lung at 5...

Note

For research use only

Application notes

BMPER antibody can be used for detection of BMPER by ELISA at 1:12500. BMPER antibody can be used for detection of BMPER by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

76 kDa

Uniprot ID
[Q8N8U9](#)
NCBI
[NP_597725](#)
Dilution Range

BMPER antibody can be used for detection of BMPER by ELISA at 1:12500. BMPER antibody can be used for detection of BMPER by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

