

Product Datasheet

UBE2F Antibody (orb1240777)


Description	UBE2F Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human UBE2F.
Target	UBE2F
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store UBE2F antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	UBE2F antibody can be used for detection of UBE2F by ELISA at 1:312500. UBE2F antibody can be used for detection of UBE2F by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	Q969M7
NCBI	NP_542409
Dilution Range	UBE2F antibody can be used for detection of UBE2F by ELISA at 1:312500. UBE2F antibody can be used for detection of UBE2F by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Mouse Brain at 2 ...



Antibody used in WB on Human HeLa at 0.2...



Antibody used in IHC on Human Brain at: ...