


Product Datasheet

CBLN4 Antibody (orb1240780)

Description	CBLN4 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	500
Conjugation	Unconjugated	400
Tested Applications	ELISA, WB	300
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CBLN4.	200
Target	CBLN4	Antibody used in WB on Human HeLa at 1 u...
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	Antibody used in IHC on Human cortex at ...
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CBLN4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in IHC on Human cortex.
Note	For research use only	
Application notes	CBLN4 antibody can be used for detection of CBLN4 by ELISA at 1:312500. CBLN4 antibody can be used for detection of CBLN4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	22 kDa	
Uniprot ID	Q9NTU7	
NCBI	NP_542184	
Dilution Range	CBLN4 antibody can be used for detection of CBLN4 by ELISA at 1:312500. CBLN4 antibody can be used for detection of CBLN4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	