



Product Datasheet

MBNL2 Antibody (orb1240866)

Description	MBNL2 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Brain at 0...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MBNL2.	
Target	MBNL2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MBNL2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Human Hela at 1 u...
Note	For research use only	
Application notes	MBNL2 antibody can be used for detection of MBNL2 by ELISA at 1:1562500. MBNL2 antibody can be used for detection of MBNL2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	28 kDa	
Uniprot ID	O95205	Antibody used in WB on Hum. Fetal Brain ...
NCBI	NP_005748	
Dilution Range	MBNL2 antibody can be used for detection of MBNL2 by ELISA at 1:1562500. MBNL2 antibody can be used for detection of MBNL2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	