

Product Datasheet

PHF6 Antibody (orb1240885)

Description

PHF6 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PHF6.

Target

PHF6

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PHF6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PHF6 antibody can be used for detection of PHF6 by ELISA at 1:312500. PHF6 antibody can be used for detection of PHF6 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

41 kDa, 35 kDa

Uniprot ID

[Q8IWS0](#)

NCBI


[NP_115834](#)

Dilution Range

PHF6 antibody can be used for detection of PHF6 by ELISA at 1:312500. PHF6 antibody can be used for detection of PHF6 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human Jurkat 5.0 ...



Antibody used in WB on Human, Mouse at 1...