

Product Datasheet

TIPIN Antibody (orb1240953)

Description	TIPIN Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Jurkat 5.0 ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TIPIN.	
Target	TIPIN	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TIPIN antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	TIPIN antibody can be used for detection of TIPIN by ELISA at 1:1562500. TIPIN antibody can be used for detection of TIPIN by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	34 kDa	
Uniprot ID	Q9BVW5	
NCBI	NP_060328	
Dilution Range	TIPIN antibody can be used for detection of TIPIN by ELISA at 1:1562500. TIPIN antibody can be used for detection of TIPIN by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	