

Product Datasheet

CASP8AP2 Antibody (orb1240956)

Description	CASP8AP2 Antibody	 <p>Antibody used in WB on Human HeLa at 0.2...</p>
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CASP8AP2.	
Target	CASP8AP2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CASP8AP2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	CASP8AP2 antibody can be used for detection of CASP8AP2 by ELISA at 1:12500. CASP8AP2 antibody can be used for detection of CASP8AP2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	223 kDa	
Uniprot ID	Q9UKL3	
NCBI	NP_036247	
Dilution Range	CASP8AP2 antibody can be used for detection of CASP8AP2 by ELISA at 1:12500. CASP8AP2 antibody can be used for detection of CASP8AP2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	