

Product Datasheet

FCRL1 Antibody (orb1240981)

Description

FCRL1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HepG2 at
2....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human FCRL1.

Target

FCRL1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store FCRL1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

FCRL1 antibody can be used for detection of FCRL1 by ELISA at 1:1562500. FCRL1 antibody can be used for detection of FCRL1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

47 kDa

Uniprot ID

[Q96LA6](#)

NCBI

[NP_443170](#)

Dilution Range

FCRL1 antibody can be used for detection of FCRL1 by ELISA at 1:1562500. FCRL1 antibody can be used for detection of FCRL1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.