

Product Datasheet

JPH1 Antibody (orb1241025)

Description	JPH1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human JPH1.	Antibody used in WB on Human Jurkat 5.0 ...
Target	JPH1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store JPH1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	JPH1 antibody can be used for detection of JPH1 by ELISA at 1:1562500. JPH1 antibody can be used for detection of JPH1 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	72 kDa	
Uniprot ID	Q9HDC5	
NCBI	NP_065698	
Dilution Range	JPH1 antibody can be used for detection of JPH1 by ELISA at 1:1562500. JPH1 antibody can be used for detection of JPH1 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	