
Product Datasheet

RABGGTA Antibody (orb1241037)

Description

RABGGTA Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody used in WB on Human Jurkat 0.2-...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RABGGTA.

Target

RABGGTA

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store RABGGTA antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in WB on Human Hela at 1 u...

Note

For research use only

Application notes

RABGGTA antibody can be used for detection of RABGGTA by ELISA at 1:62500. RABGGTA antibody can be used for detection of RABGGTA by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

65 kDa

Uniprot ID
[Q92696](#)
NCBI
[NP_878256](#)
Dilution Range

RABGGTA antibody can be used for detection of RABGGTA by ELISA at 1:62500. RABGGTA antibody can be used for detection of RABGGTA by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

 100kDa
80kDa

 60kDa
40kDa
20kDa

RABGGTA

 80
60
40
20

 Rabbit Anti-RABGGTA
Sample Type: Human Hela
Antibody Concentration: 1µg/ml