

Product Datasheet

TRMT2B Antibody (orb1241070)

Description	TRMT2B Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human HepG2 at 2....
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CXorf34.	
Target	TRMT2B	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CXorf34 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	CXorf34 antibody can be used for detection of CXorf34 by ELISA at 1:1562500. CXorf34 antibody can be used for detection of CXorf34 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	56 kDa	
Uniprot ID	Q96GJ1	
NCBI	NP_079193	
Dilution Range	CXorf34 antibody can be used for detection of CXorf34 by ELISA at 1:1562500. CXorf34 antibody can be used for detection of CXorf34 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	