


Product Datasheet

PAPD7 Antibody (orb1241096)

Description	PAPD7 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human HepG2 at 0....
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human POLS.	
Target	PAPD7	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store POLS antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	POLS antibody can be used for detection of POLS by ELISA at 1:312500. POLS antibody can be used for detection of POLS by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	60 kDa	
Uniprot ID	Q5XG87	
NCBI	NP_008930	
Dilution Range	POLS antibody can be used for detection of POLS by ELISA at 1:312500. POLS antibody can be used for detection of POLS by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	