


Product Datasheet

PFN1 Antibody (orb1241141)

Description	PFN1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PFN1.
Target	PFN1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PFN1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PFN1 antibody can be used for detection of PFN1 by ELISA at 1:12500. PFN1 antibody can be used for detection of PFN1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	15 kDa
Uniprot ID	P07737
NCBI	NP_005013
Dilution Range	PFN1 antibody can be used for detection of PFN1 by ELISA at 1:12500. PFN1 antibody can be used for detection of PFN1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.




Antibody used in WB on Human Liver at 0....



See Immunoblot 2 Data and Customer Feedback for more information.

Antibody used in WB on Human at: 1:2000 ...



See Immunoblot 3 Data and Customer Feedback for more information.

Antibody used in WB on Mouse at: 1:2000 ...