
Product Datasheet

GLYCTK Antibody (orb1241208)

Description

GLYCTK Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

 Antibody
used in WB
on Human
HepG2 at
1....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GLYCTK.

Target

GLYCTK

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GLYCTK antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GLYCTK antibody can be used for detection of GLYCTK by ELISA at 1:312500. GLYCTK antibody can be used for detection of GLYCTK by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

55 kDa

Uniprot ID
[Q8IVS8](#)
NCBI
[NP_660305](#)
Dilution Range

GLYCTK antibody can be used for detection of GLYCTK by ELISA at 1:312500. GLYCTK antibody can be used for detection of GLYCTK by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.