
Product Datasheet

SV2A Antibody (orb1241257)

Description	SV2A Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Brain at 0...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SV2A.	
Target	SV2A	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SV2A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	SV2A antibody can be used for detection of SV2A by ELISA at 1:62500. SV2A antibody can be used for detection of SV2A by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	83 kDa	
Uniprot ID	Q7L0J3	
NCBI	NP_055664	
Dilution Range	SV2A antibody can be used for detection of SV2A by ELISA at 1:62500. SV2A antibody can be used for detection of SV2A by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	