
Product Datasheet

FCGRT Antibody (orb1241323)

Description	FCGRT Antibody	
Species/Host	Rabbit	70kDa
Reactivity	Human	40kDa
Conjugation	Unconjugated	30kDa
Tested Applications	ELISA, WB	20kDa
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human FCGRT.	Antibody used in WB on Human HepG2 at 0....
Target	FCGRT	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and may contain up to 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store FCGRT antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	FCGRT antibody can be used for detection of FCGRT by ELISA at 1:62500. FCGRT antibody can be used for detection of FCGRT by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	40 kDa	
Uniprot ID	P55899	
NCBI	NP_004098	
Dilution Range	FCGRT antibody can be used for detection of FCGRT by ELISA at 1:62500. FCGRT antibody can be used for detection of FCGRT by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	