

Product Datasheet

CCT8 Antibody (orb1241396)


Description	CCT8 Antibody
Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CCT8.
Target	CCT8
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CCT8 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	CCT8 antibody can be used for detection of CCT8 by ELISA at 1:62500. CCT8 antibody can be used for detection of CCT8 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	Q5RAP1
NCBI	NP_006576
Dilution Range	CCT8 antibody can be used for detection of CCT8 by ELISA at 1:62500. CCT8 antibody can be used for detection of CCT8 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Jurkat 0.2-...



Antibody used in WB on Human, Mouse.



Antibody used in WB on Human HEK293T cel...