
Product Datasheet

CHAF1B Antibody (orb1241397)

Description

CHAF1B Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CHAF1B.

Target

CHAF1B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CHAF1B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CHAF1B antibody can be used for detection of CHAF1B by ELISA at 1:62500. CHAF1B antibody can be used for detection of CHAF1B by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW


61 kDa

Uniprot ID
[Q13112](#)
NCBI
[NP_005432](#)
Dilution Range

CHAF1B antibody can be used for detection of CHAF1B by ELISA at 1:62500. CHAF1B antibody can be used for detection of CHAF1B by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 5.0 ...